



HAVE YOU CHECKED THE FOLLOWING BEFORE SENDING OUT A SAMPLE?

Our Specimen Rejection Policy

Specimens would be rejected in the following cases:




- Unlabelled specimen




- Specimen sent in wrong container
- Inappropriate container
- Wrong preservative for test



- Incorrectly labelled specimens
- (Mismatch Specimen and Form)



- Insufficient sample volume
(Volume not to the mark)
e.g Citrated blood



- Broken
- Leaking container
- Spilt sample



- Contaminated sample or Contaminated request form



- Specimen received in syringe

Help us, help you to avoid delays or rejected specimens

International CAP Accreditation 2013 • Bringing innovative new tests to Singapore • Certified Laboratory Excellence

innovative
DIAGNOSTICS
leading the way

cap
ACCREDITED™